

Automatization – Macro Recording

Introductory Course to Image Analysis in Life Science



LECLERC SIMON

Life Line

- <https://imagej.net/ij/>
- <https://imagej.net/ij/developer/macro/macros.html>
- <https://imagej.net/ij/developer/macro/functions.html>


Material for the exercises:

Making use of Macros to automatize image processing



The Basics


Open the Macros Editor: **Plugin > New > Macro**



The Basics: Variables

All variables needs to be declared

A lot of mathematical functions are there




The screenshot shows a Java IDE window titled "*Macro.iqm (Running)". The code in the editor is as follows:

```
1 // Initialization
2 a = 6;
3 b = 3.0;
4 c = "one";
5
6 // some operation
7 d = a + b;
8
9 // result
10 print(d);
```

The line `print(d);` is highlighted in yellow. To the right of the code editor is a "Log" window with a menu bar (File, Edit, Font) and the following output:

```
Hello World
9
```



The screenshot shows a Java IDE window titled "*Macro.iqm". The code in the editor is as follows:


```
1 // Initialization
2 a = 6;
3 b = 3.0;
4
5 // Mathematical print
6 print(pow(a,b)); //a to the power of b
7 print(sqrt(a)); // square root
8 print(round(b)); // round
9 print(abs(-b)); // absolute
```

The line `print(abs(-b)); // absolute` is highlighted in yellow. To the right of the code editor is a "Log" window with a menu bar (File, Edit, Font) and the following output:

```
216
2.4495
3
3
```

The Basics: String Manipulation

Concatenate, select a substring, replace a part of the string, check if a string starts or ends with a selection



The screenshot shows a Java IDE with a file named 'Macro.ijsm'. The code on the left performs various string operations: initialization of 'a' and 'b', concatenation, substring extraction, replacement, and checks for start/end characters. The output on the right shows the results of these operations: 'Hello World', a blank line, 'Worm', and two occurrences of the number '1'.

```
Macro.ijsm
1 // Initialization
2 a = "Hello";
3 b = "World";
4
5 // string manipulation print
6 print(a+" "+b); // concatenate
7 print(substring(a, 2, 4)); // substring
8 print(replace(b, "ld", "m")); // replace
9 print(startsWith(a, "He")); // check the start of the string
10 print(endsWith(b, "ld")); // check the end of the string
```

File Edit Font

Hello World


Worm

1

1

Autocompletion

Use the Autocompletion to simplify your life!



The Recorder

Easiest way to do a Macro: Record the actions that you do in Fiji!

Open it with: `Plugins > Macros > Record...`

- Open the Blobs image (`File > Open Sample > Blobs`)
- Run a Gaussian Blur with sigma of 3 (`Process > Gaussian Blur...`)
- Set a Otsu threshold (`Image > Adjust > Threshold...`)
- Set the parameters to measure as area, centroid and shape (`Analyze > Set Measurements...`)
- Analyze the segmented image (`Analyze > Analyze Particle...`)
- Then on the Recorder, click on Create to automatically generate the macro

The Recorder

```
1 run("Blobs (25K)");
2 run("Gaussian Blur...", "sigma=3");
3 setAutoThreshold("Otsu");
4 //run("Threshold...");
5 //setThreshold(123, 255);
6 setOption("BlackBackground", true);
7 run("Convert to Mask");
8 run("Close");
9 run("Set Measurements...", "area centroid shape redirect=None decimal=3");
10 run("Analyze Particles...", "size=0-Infinity display clear");
11
```

Line starting by ‘//’
are comments,
piece of text that will
NOT be executed

The Recorder

```
1 //open the image
2 run("Blobs (25K)");
3
4 //processing
5 run("Gaussian Blur...", "sigma=3");
6
7 //Segmentation
8 setAutoThreshold("Otsu");
9 setOption("BlackBackground", true);
10 run("Convert to Mask");
11
12 //Analyze
13 run("Set Measurements...", "area centroid shape redirect=None decimal=3");
14 run("Analyze Particles...", "size=0-Infinity display clear");
15
```

The 'run' function

- You may have observed a pattern in the first macro: the run function is called multiple times. And all the parameter(s) of this function is string.
- The first parameter is the name of the process called, such as Convert to Mask.
- In the case the name ends with ..., this means that this function require some additional parameters, which are describe in the second string, such as sigma=3 for the Gaussian Blur... by example.
- In this case, string manipulation, and mainly concatenation, is very important to be able to dynamically modify some parameters.

```
run("Convert to Mask");
```

```
run("Gaussian Blur...", "sigma=3");
```

Exercise 1: Make a Simple Montage

- Open an Image (2_22.tif for example)
- Record the following actions:
 - Do a Z Projection, Max Intensity (Image > Stacks > Z-Project...)
 - Rename the window (right click, Rename...)
 - Adjust Contract & Brightness (Image > Adjust > Brightness/Contrast...)
 - Give a colour (Image > Lookup Table > Cyan)
 - Repeat the last 2 points after changing Colour
 - Make a Montage (Image > Stacks > Make Montage...)
 - Create the Macro


Does EVERY actions get recorded? Try your Macro on another Image!

Colour: Composite Image

You can make composite image with the **Image > Color > Channels Tools...**

Select 'Composite' and unselect Channel 3.

Then **More, Create RGB** to get flatten image, then you can add a scale bar
Analysis > Tools > Scale Bar...



Window Selection

Image operation always happen on the **currently selected image**.

It is possible to select window with the 'selectImage("image_name")' function.

Some Image operation will create a new Image (Z Projection).

Other will apply the operation in place (Gaussian Blur).


In case of doubt, duplicate the image (Right click, Duplicate...)

Image Type

- Image type is indicated at the top left of the image

Some operations can only be made on some type!
8 and 16 bit are integer
32 bit are float
RGB are for display and saving

Copy/Paste Slice is possible between identical image type



Exercise 2: Make a Better Montage

- Open an Image (2_22.tif for example)
- Starting from the previous macro result
- Record the following actions:
 - With the Channel tools, make a composite image with only the fluorescent channels
 - Convert the image as RGB, add a scale bar, and copy it **Edit > Copy**
 - Select the Zprojection
 - Add a slice to it **Image > Stacks > Add Slice** and paste the composite **Edit > Paste**
 - Make the Montage
 - Save as PNG **File > Save As > PNG**

Variables: Getting Name

- To keep track of your images in the macro, it is important to rename, but also getting the name of the currently selected image

`getTitle();` return the name of the currently open image

`Img_name = getTitle();`

It is possible to get much more than the title if required. Try the 'get' and check the autocompletion!

Exercise 3: Keep montage!

- Open an Image (2_22.tif for example)
- Starting from the previous macro result
- Modify the macro to:
 - Get the title of the initial image
 - Print it in the log window
 - When saving the montage as PNG, adjust the original name

Run your macro on multiple image and be sure that the montage image name is unique (not erased)


Programming: For Loop

- This can be translated as 'Do this job X times', or for example 'Stir the pot 10 times'.
- In the ImageJ Macro language, it requires a little more information:
 - An initialization stage (where we start, at 0 stir)
 - A condition (when to stop, at 10 stir)
 - An increment/decrement (how to change each time, 1 stir at a time)
 - The code to be executed in each iteration (stir)

Programming: For Loop

```
for (i = 0; i < 10; i++) {  
    print(i)  
}
```

If the condition is always false = infinite loop!



Exercise 4: Running the Macro through all images

- Starting from the previous macro result
- Modify the macro to:
 - Loop through all files within a directory – check ‘for’ and autocompletion!
 - You will need to assign a directory variable
 - Include the previously written code to make the montage
 - Open the image using the Bio-Format Importer **Plugin > Bio-Format**

Test the macro and note the issue(s)

Programming: if condition

- Check if the condition is true

```
if (a == b) {  
    print(a+' is equal to '+b);  
}  
else {  
    print(a+' is NOT equal to '+b);  
}
```

Programming: if condition

It is possible to chain condition check with the keyword else if

It is possible to check multiple condition at the same time with
&& (and), and || (or)

Some function return a boolean (true/false) by default:

is(« grayscale »);

Exercise 5: Running the Macro through tif images

- Starting from the previous macro result
- Modify the macro to:
 - Check if the file that we will open is a tif
 - If it is, process as we have been doing
 - Else go to the next image

Programming: testing

It is important to test that the macro is working on multiple but similar images

```
run("Close All");
```

```
print("\\Clear");
```

```
run("Clear Results");
```

To start a macro with a clean slate

Programming: Array

The macro language can store data in Array

```
color_array = newArray("Cyan", "Yellow", "Grays", "Green",  
"Magenta");
```

It is possible to loop through an Array with a for loop, or modify and specific entry by knowing its index

```
color_array[0] = "Blue"
```

Check the 'Array' autocompletion for more information

Programming: GUI

The macro language also allow you to build a full GUI with interactive window.

But before that, there is 2 interesting functions:

`getDirectory("Choose directory");` // Open a GUI window to prompt the user to select a directory

`waitForUser(string);` // Stop the macro and wait for the user to click on 'OK' before proceeding.

Programming: Batch processing

Displaying images has a computational cost, which can slow down or even freeze Window!

It is possible to hide all images to avoid this extra cost by using:

- **setBatchMode("hide")**: Hide active image
- **setBatchMode("show")**: Display active image
- **setBatchMode("exit and display")**: Display ALL images

Optional Exercise: More general Montage

- Starting from the previous macro result
- Modify the macro so that:
 - The Lookup Table section is more polyvalent
 - `getDimensions`
 - `Stack.setChannel` within a for loop
 - Adjust inactive channel for Grays
 - Add a 'choose directory' window and BatchMode
 - Super optional: Add a GUI to allow identification and colour selection of the brightfield channel



UNIVERSITY OF
GOTHENBURG

Thanks for your attention